# **PRESIDENT System 360**

## **heavy body • surface activated**

2	Gebrauchsinformation
5	Instructions for use
8	Mode d'emploi
11	Información para el uso
14	Informazioni per l'uso
17	Bruksanvisning
20	Gebruiksaanwijzing
23	Brugsanvisning
26	Οδηγίες χρήσεως
29	Instruções de Utilização
32	Kullanma Talimatı

**coltène**   
**whaledent**<sup>®</sup>

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### Definition

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## Instructions for use (English)

### Definition

Coltène® President System 360 heavy body is a dental impression material for use in dentistry with dynamic mixing tips for usage in commercial automatic mixing machines.

### Material type

Polyvinylsiloxane, addition-type, surface activated silicone elastomer, medium and high viscosity.

### President heavy body

Base: light brown  
Catalyst: caramel  
– ISO 4823, Type 1,  
high consistency

### Indications President heavy body

- Tray material for the 1-step or double mix simultaneous technique
- Tray material for primary impression in 2-step «heavy body/wash» technique (similar to putty/wash technique)
- Tray material for the dual arch bite technique (closed mouth triple tray technique)
- For full and partial denture impressions
- For all crowns and bridges, border moldings, transfer copings or «pick-ups» and implant cases (once healing has taken place and case is ready for impression taking)

### Clinical time heavy body

### Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

### Side effects / Interactions

Polyvinylsiloxanes have very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient. However, they should remain in the mouth not more than twice the recommended setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

### Trays

Tray selection depends on impression technique and professional preference of metal or rigid plastic stock trays, double arch bite trays or custom trays. For perfect adhesion, we recommend applying a thin layer of Coltène® Adhesive (1 min drying) or any other brand of adhesive specified for use with polyvinylsiloxane impression materials. Do not use adhesive for polyether materials. Apply the proper adhesive to the metal or plastic tray. Do not apply adhesive to the gauze web section of double arch bite trays!

### Preparation

1. Remove the white safety bolt from the large foil bag by turning it clockwise.  
**Attention: Only use base and catalyst with identical **LOT** number.**
2. Set up the supporting cartridge. Release the safety clip.
3. Push both foil bags with the premounted adapter heads into the supporting cartridge. Take care to ensure that the cartridge body and the

adapter head notches are aligned. This ensures optimum alignment of the outlet channels.

4. Press the adapter heads with the ball of the thumb onto the supporting cartridge until they lock firmly into final position. The foil bags are automatically pierced by the pin pressing the adapter heads onto the cartridge.
5. Insert the filled cartridge body into the mixing machine.
6. Switch on the device and extrude sufficient material from the newly pierced tubes to obtain even layer.
7. Assemble the dynamic mixing tip. Take care to ensure the right alignment and the engagement of the drive shaft. Reposition the safety clip. Fill the tray. The dynamic mixing tip should always stay in the material to avoid air bubbles.
8. Immediately after the tray filling procedure remove the dynamic mixing tip and ease the pressure of the automatic mixing machine.
9. Put a new dynamic mixing tip and do not remove it from the cartridge as it acts as a seal.
10. Before using again, check that the outlets of the activating heads are not blocked. If necessary clean the outlets with a pointed instrument and check for proper extrusion of catalyst and base every time and continue as usual.
11. Remove the empty foil bags including the adapter heads from the supporting cartridge and dispose of them.

**Important:** In order to ensure optimum performance, please use only Coltène Whaledent® dynamic mixing tips.

**Recommendations for use «double mix technique»**

First fill the tray with Coltène® President System 360 heavy body. As soon as tray is halfway filled, begin simultaneously to syringe material around preparation using either a hand held impression syringe or the Coltène® President light or regular body material. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place tray into the patient's mouth. Press for 2–3 s and hold in position.*

**Check intraorally if the material is set before removing from the mouth.**

**Please note the following:** If a 2-step «heavy body/wash» technique is used, the primary impression made with heavy body material must be carefully cleaned and dried before use. This will help to guarantee good adhesion between the wash and tray material.

Some latex gloves and surfaces contaminated by these gloves (teeth, preparations, retraction cords, etc.) may interfere with the setting process of polyvinylsiloxanes. It may be necessary to use vinyl gloves or to thoroughly wash, rinse and dry hands (wash for 15 s with a detergent, rinse with water for another 15 s and dry). Compounds containing eugenol or hemostatic substances can also impede perfect setting. If hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation.

**Disinfection**

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

**Fabrication of models**

The stone models can be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (tested for 7 days). Coltène® impression materials are compatible with all dental stones (e.g. Fuji-rock, Vel-mix), epoxy resin and polyurethane.

**Electro plating**

Compatible with all commercially available copper plating and silver plating baths.

**Cleaning of trays**

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or petroleum ether will dissolve the Coltène® Adhesive. These solvents should only be used in a well-ventilated area. The trays can also be ultrasonically cleaned and sterilized.

**Shelf life and storage**

Recommended storage of President System 360 heavy body material is 15–23 °C / 59–73 °F and 50% relative air humidity. Avoid exposure to direct heat and sunlight. Coltène® President System 360 heavy body fulfill the intended purpose at least until the expiry date displayed on the package and cartridge.

**Caution**

Federal law restricts this device to sale by or on the order of a dentist.

**Marking**

The expiry date and **LOTI** number are shown on the package and cartridge.

**Technical data ISO 4823:2000**

Measurements are based on 23 °C / 73 °F room temperature and 50% relative humidity.

**President System 360 heavy body**

<u>Mixing time:</u>	15 s
<u>Total working time:</u>	45 s
<u>Oral setting time:</u>	180 s
<u>Elastic recovery:</u>	99.75 %
<u>Strain in compression:</u>	2.3 %
<u>Dimensional change:</u>	- 0.18 %

**Date of issue**

April 2007

Only supplied to dentists and dental laboratories or upon their instructions.

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### President heavy body

Base: brun clair  
Catalyseur: caramel  
– ISO 4823, Type 1,  
high consistency

### Indications de President heavy body

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**Coltène/Whaledent AG (Manufacturer)**