coltène®

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Speedex light body • surface activated

Definition

Speedex light body is a silicone-based impression material for use in dentistry.

Material type

Polysiloxane, condensation-type silicone elastomer.
 Low viscosity
 Base: dark blue
 Universal Activator (separate): green
 – ISO 4823, Type 3, low consistency

Indications

- Wash material for the putty-wash impression technique
- Syringe material for the simultaneous mix technique
- Impression material for relining

Clinical time (normal dosage)

≥ 4:30 min		
≤ 1:30 min	≤ 1:00 min	≥ 3:00 min
		
≤ 0:30 min	≤ 1:00 min	≥ 3:00 min
23 °C / 73 °F		

Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions


Polysiloxanes have very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

In case of contact of the catalyst with the eyes, wash out with plenty of water and then consult an ophthalmologist.


Trays

Rigid stock or individual trays. For perfect adhesion, we recommend to apply a thin layer of Coltène® Adhesive to all trays.

Dosage


Squeeze out equal strand lengths of base material and activator.

Always dose universal activator first for better results.
Close tube immediately after use!

Over- and underdosage:

By over- and underdosing the universal activator setting time can be accelerated or slowed down.


Clinical time

(Overdosage, 20% more universal activator)

≥ 3:00 min		
≤ 1:00 min	≤ 0:40 min	≥ 2:00 min
		
≤ 0:20 min	≤ 0:40 min	≥ 2:00 min
23 °C / 73 °F		

Even if you chose the shortest setting time, there will be enough working time to apply the material correctly and seat the tray.

Mixing

Take up one component with spatula and spread evenly on the other component.

Impression

Insert the tray and press it into position (2–3 s). Hold it without pressure until material is set. Material residue needs considerably more time to set at room temperature. You should therefore check intraorally if the material is set before removing from the mouth. Intensive mixing, high temperatures or **overdosage of the universal activator speed up** the setting; low temperatures and **underdosage of the universal activator slow down** the setting.

Important: In order to guarantee good cohesion with the wash material, the primary impression must be carefully cleaned and dried before further use. If H₂O₂ hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation.

Disinfection

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

Manufacture of models

The ideal time is between 30 min and 72 h after taking the impression. Thanks to its **extraordinary dimensional stability, models can be poured up to 7 days after impression taking without any detrimental effects.** The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. All industry-standard dental stone model materials (i.e. Fujirock Dental Stone, Hard Rock Dental Stone) can be used.

Electro plating

Compatible with all commercially available copper plating and silver plating baths.

Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or light petrol will dissolve the Coltène® Adhesive. These solvents should only be used in a well-ventilated area. The trays can be cleaned and disinfected/sterilized as usual.

Shelf life and storage

Speedex light body fulfills the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed tubes at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

Marking

The expiry date and **LOT** number are shown on the package and cartridge.

Technical data ISO 4823:2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

Mixing time:

0:30 min

Oral setting time:

3:00 min

Instructions for use date

February 2009

Only supplied to dentists and dental laboratories or upon their instructions.

Caution: Federal law restricts this device to sale by or on the order of a dentist.

Coltène/Whaledent AG (Manufacturer)

Feldwiesenstrasse 20
 9450 Altstätten/Switzerland
 Tel. +41 (0)71 757 53 00
 Fax +41 (0)71 757 53 01
 info@coltenewhaledent.ch

For MSDS see

www.coltenewhaledent.com

Speedex light body • surface activated

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