coltène® Speedex putty

Definition

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coltène® Speedex putty

Definition

Speedex putty is a silicone-based impression material for use in dentistry.

Material type

Polysiloxane, condensation-type silicone elastomer, high viscosity putty material.

Base: light grey


Universal activator (separate): green

– ISO 4823, Type 1, high consistency

Indications

- Primary impression in the putty-wash impression technique
- Tray material for the simultaneous mixing technique
- Tray material for the copperband impression technique
- Impressions for study models, orthodontic models, matrices, protective insulation of teeth during flasking of the denture

Setting time – normal dosage

Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions


Polysiloxanes have very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

In case of contact of the catalyst with the eyes, wash out with plenty of water and then consult the ophthalmologist.

Trays

Rigid stock trays. For perfect adhesion, we recommend to apply a thin layer of Coltène® Adhesive to all trays.

Dosage

Dose Speedex putty with measuring-scoop (scrape off surplus). Spread out the dosed amount in the palm of hand.

For each scoop base press the scoop-rim once into the mass. Dose one strand-length of universal activator per circle-diameter.


Close tube of universal activator immediately after use!

Over- and under dosage:

By over- and underdosing the universal activator setting time can be accelerated or slowed down.

Setting time – Overdosage 20% more universal activator

Even if you chose the shortest setting time, there will be enough working time to apply the material correctly and seat the tray.

Mixing

Fold material with fingertips and knead vigorously until an absolutely uniform color is attained, but for at least 30 s. When using mechanical mixing devices, observe the instructions of the manufacturer.

Impression

Insert the tray and press it into position (2–3 s). Hold it without pressure until material is set. Material residue needs considerably more time to set at room temperature. You should therefore check intraorally if the material is set before removing from the mouth. Intensive kneading, high temperatures or **overdosage of the universal activator speed up** the setting; low temperatures and **underdosage of the universal activator slow down** the setting.

Important: In order to guarantee good cohesion with the wash material, the primary impression must be carefully cleaned and dried before further use.

If H₂O₂ hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with luke-warm water afterwards in order to avoid bubble formation.

Disinfection

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

Manufacture of models

The ideal time is between 30 min and 72 hours after taking the impression. The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in luke-warm clear water afterwards. All industry-standard dental stone model materials (i.e. Fujirock Dental Stone, Hard Rock Dental Stone) can be used.

Electro plating

Compatible with all commercially available copper plating and silver plating baths.

Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or light petrol will dissolve the Coltène® Adhesive. These solvents should only be used in a well ventilated area. The trays can be cleaned and disinfected as usual.

Shelf life and storage

Speedex putty fulfills the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

Marking

The expiry date and **LOT** number are shown on the package and cartridge.

Technical data ISO 4823:2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

Mixing time:

0:30 min

Oral setting time:

3:00 min

Instructions for use date

February 2009

Only supplied to dentists and dental laboratories or upon their instructions.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

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Coltene/Whaledent AG (Manufacturer)