DENTSPLY

Aquasil Ultra XLV or Aquasil Ultra LV Regular Set

**Smart Wetting® Impression Material in cartridge
or digit™ Targeted Delivery System**

CE

TABLE OF CONTENTS

English	2
Español	17
Français	35
Deutsch	53
Italiano	71
Português	89
Dansk	107
Nederlands	123
Svenska	141
Ελληνικά	157

DIRECTIONS FOR USE – ENGLISH

CLASSIFICATION

Quadrafunctional Hydrophilic Addition Reaction Silicone Impression Material
Type 3: Light-Bodied Consistency

Caution: U.S. federal law restricts this device to sale by or on the order of a dentist. For dental use only.

DESCRIPTION

Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Smart Wetting® Impression Material is a quadrafunctional hydrophilic addition reaction silicone, light body, elastomeric impression material with excellent hydrophilic properties, dimensional accuracy, high tear strength, and resistance to permanent deformation. Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material with Smart Wetting® improves both wetting on the tooth and tissue surface and model detail reproduction. Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material is available in 50mL cartridge delivery and unit dose digit™ Targeted Delivery System cartridges. It complies with the requirements of ISO 4823 for dental elastomeric impression materials.

Aquasil Ultra Impression Material is to be used following the application of a suitable tray adhesive such as Caulk® Tray Adhesive or Silfix Tray Adhesive. Use of other tray adhesives with Aquasil Ultra Impression Material is at the discretion and sole responsibility of the dental practitioner.

COMPOSITION

Polydimethyldivinylsiloxane polymer; Silicon Dioxide; Polymethylhydrogen Siloxane; Surfactant; Plasticizer; Pigments

INDICATIONS FOR USE

Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material is suitable for: All tooth, implant (closed flap) or tissue impression techniques where a light body (wash) material would be desired by the operator.

1. Dual phase, one-step full arch impression technique with Aquasil Ultra Monophase Regular Set, and Aquasil Ultra Heavy Regular Set (cartridge or DECA™ dispensed), Aquasil Putty, or Reprosil® Regular, Heavy or Putty Impression Materials.
2. Dual phase, one-step quadrant arch impression technique with Aquasil Ultra Rigid Regular Set, Aquasil Ultra Monophase Regular Set, Aquasil Ultra Heavy Regular Set (cartridge or DECA™ dispensed), Aquasil Putty, or Reprosil® Regular, Heavy or Putty Impression Materials.
3. Dual phase, one step double arch (Caulk First Bite® Trays) impression technique with Aquasil Ultra Rigid Regular Set, Aquasil Ultra Monophase Regular Set or Aquasil Ultra Heavy Regular Set (cartridge or DECA™ dispensed) or Reprosil® Regular or Heavy Impression Materials.
4. Dual phase, two-step full or partial arch impression technique with Aquasil Putty or Reprosil® Putty Impression Materials.
5. Edentulous and abutment/tissue overdenture impression with properly fitted tray.
6. Dual phase, any combination, for the duplication of models.

Classification ISO 4823	Aquasil Ultra XLV: Type 3: Light-bodied consistency Aquasil Ultra LV: Type 3: Light-bodied consistency Aquasil Ultra Cartridge & DECA™ Monophase: Type 2: Medium-bodied consistency Aquasil Ultra Cartridge & DECA™ Heavy: Type 2: Medium-bodied consistency Aquasil Ultra Rigid: Type 1: Heavy-bodied consistency
Working Time*(22°C/72°F)	2'15-2'45" (room temperature mixing and preparation)
Working Time*(37°C/98°F)	Maximum 1'10" (intraoral syringing prior to tray insertion)
% Recovery from deformation	>98.0%
% Strain in compression	Aquasil Ultra XLV: 4-5 Aquasil Ultra LV: 3-4 Aquasil Ultra Cartridge & DECA™ Monophase: 3-4 Aquasil Ultra Cartridge & DECA™ Heavy: 2-3 Aquasil Ultra Rigid: 1-2
Detail Reproduction	< 20 microns
Linear dimension change	< 0.50%
Additional information:	
Maximum intraoral Working Time (tray insertion time)	1'10" from start of syringing
Setting Time (Mouth Removal Time)	5'00" from start of mix
Mixing ratio by volume	1 part base to 1 part catalyst

*Ta instruments AR-1000 Oscillating Rheometer
 Data on File

CONTRAINDICATIONS

None Known.

WARNINGS

1. Avoid prolonged or repeated exposure of the Aquasil Ultra impression materials with skin and eyes. Irritation and possible corneal damage may result. Skin rash, oral mucosa irritation, or other allergic reactions may result in susceptible individuals.
Eye and Skin Contact: Flush eyes with flowing water for 15 minutes and consult a physician. Flush skin with flowing water for 15 minutes, then wash area with soap and water after contact. Flush oral tissues with copious amounts of water. Consult a physician if rash persists.
Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.
2. Do not use Aquasil Ultra Impression Material as a temporary reliner. Use of the material in this manner may cause irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.

PRECAUTIONS

1. This product is intended to be used only as specifically outlined in the *Directions For Use*. Any use of this product inconsistent with the *Directions For Use* is at the discretion and sole responsibility of the practitioner.
2. Aquasil Ultra Impression Material should not be intermixed with, or used in conjunction with other manufacturer's vinyl polysiloxanes. Aquasil Ultra Material should not be intermixed with, or used in conjunction with polyether, polysulfide or conventional (condensation cured) silicones.
3. Aquasil Ultra Regular Set Impression Materials should not be intermixed with Aquasil Ultra Fast Set Impression Materials.

4. Insufficient data exist to support the use of Aquasil Ultra Impression Materials for impressions on exposed bone. Aquasil Ultra material should only be used on or around intact mucosa.
5. Allow Aquasil Ultra Impression Material to reach room temperature prior to use. Higher temperatures reduce work times and laboratory bench set times (faster), lower temperatures increase them (slower). Clinically, intraoral set time (Mouth Removal Time) is unaffected within the range of allowable storage and initial mixing temperatures.
6. To obtain optimum physical properties, the mix should be completely homogeneous (streak free).
7. Do not contaminate retraction cords with latex or “sulfur-containing” polymeric gloves. Handle cords with college pliers. Do not touch prepared tooth with gloves.
8. Use of some hand soaps and lotions can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during mixing will avoid this possibility.
9. If the preparation is treated with hydrogen peroxide, rinse with copious amounts of water before making the impression.
10. Aquasil Ultra Impression Materials should be used with suitable impression trays pre-treated with tray adhesive (*see Step-by-Step Instructions*). Refer to tray adhesive manufacturer’s *Direction for Use* for compatibility.
11. Material should extrude easily. **DO NOT USE EXCESSIVE FORCE**. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture.
12. digit™ targeted delivery system unit dose cartridges are designed for one-time use. Discard after use. Do not reuse.
13. Due to the extreme hydrophilicity of Aquasil Ultra impression materials, use of this material in the dual phase, two-step crown and bridge impression technique is only recommended in conjunction with Aquasil Putty or Reprosil® impression materials (available separately) for primary tray impressions. Two-step technique is NOT recommended for double arch trays. When used as a primary crown and bridge impression, tray material should be covered with a polyethylene spacer sheet, and protected from direct contact with oral tissue (*see complete Directions For Use for the selected tray material*).

14. **Storage:** Store Aquasil Ultra Impression Material at or below room temperature (25°C/77°F). Refrigerated storage is acceptable when not in use. Allow material to reach room temperature prior to use. Work time is affected by temperature (see *Precautions, Step-by-Step Instructions*). Protect from moisture. Do not freeze. Do not use after expiration date.

ADVERSE REACTIONS

1. Corneal damage may result with prolonged eye exposure to the impression material. See *Warnings*
2. Allergic contact dermatitis, oral mucosa irritation, and other allergic reactions may occur in susceptible individuals. See *Warnings*
3. Skin irritation may result from direct exposure and skin defatting from prolonged exposure. See *Warnings and Precautions*
4. Bowel obstruction or other digestive distress may result from ingestion of mixed impression material. See *Warnings*

STEP-BY-STEP INSTRUCTIONS

CARTRIDGE DISPENSING GUN INSTRUCTIONS

1. Raise the release lever vertically upward, while simultaneously pulling the plunger all the way back in the dispenser handle.
2. Cartridge Loading
 - 2.1 Open cartridge lock by lifting up top clasp.
 - 2.2 Orient and insert cartridge with v-shaped notch facing down.
 - 2.3 Close the top clasp to lock cartridge into dispenser gun.
 - 2.4 Remove the cartridge cap by turning 90° counterclockwise. The cartridge cap can be replaced on cartridge when storing the cartridge after initial use or the used mixing tip can be left in place until the next use after disinfection to serve as a self-sealing cap.
 - 2.5 Dispense a small amount of base and catalyst before installing mix tip to ensure even flow from cartridge. Use gentle pressure. Be sure no plug is present. If any plug blocks flow, clear with an instrument. Wipe away excess from the

cartridge carefully so base and catalyst do not cross contaminate and cause obstruction of the nozzle.

- 2.6 Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside colored rim of the mix tip with the v-shape notch on the cartridge flange. If mix tip does not seat easily, be sure the internal end view of the mix tip is aligned correctly. The mix tip itself can turn in its cap. The two holes should line up with the v-shaped piece that is on the outside rim of the mix tip cap to facilitate inserting. Turn tip to align in its cap.
- 2.7 When mix tip colored rim v-notch is lined up with the v-shaped notch on the cartridge flange, turn colored mix tip cap 90° in a clockwise direction to lock in place on cartridge.

3. Dispenser Operation

- 3.1 Squeeze dispenser trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released.
- 3.2 The selected Aquasil Ultra syringe (wash) Impression Material may be dispensed directly to the mouth with provided mix tips. The intraoral tip can be inserted into the end of the mix tip (stator) for direct intraoral syringing

OR

An intraoral syringe may be filled by extruding the mixed syringe (wash) material directly into the syringe by back filling.

For subsequent usage, remove the used mixing tip and squeeze the dispensing gun trigger to dispense (bleed) material from the two openings of cartridge. Wipe end of cartridge. Install a new mixing tip and twist one-quarter turn to lock to position.

To remove the cartridge, hold the release lever vertically while simultaneously withdrawing plunger in a backward direction. Dispose of cartridge when empty.

DIGIT™ TARGETED DELIVERY SYSTEM UNIT DOSE DISPENSED

1. Assemble clean syringe by inserting plunger into open end of syringe barrel. Assure that plunger moves freely within syringe barrel. Have clean, assembled syringe, new mix tip and new intra-oral tip (if applicable) available. Do not assemble cartridge/mix tip until ready for use.
2. When ready for use, place the circular disk of the unit dose cartridge into the slotted end on the syringe plunger. Bend unit dose cartridge to snap off the circular disk. (Do not twist cartridge.) Disc should snap off cleanly, exposing both cartridge ports.
3. Grasp clean mix tip in one hand, and cartridge in the other. Insert unit dose cartridge into the mix tip assembly. Press firmly until the cartridge is fully seated and snapped into place. Both cartridge spurs must be fully locked into mix tip slots before proceeding. If not locked, apply additional pressure until locking, or discard mix tip and select another tip for assembly.
4. If desired for material application, attach intra-oral tip.
5. Retract syringe plunger until stop is felt. Do not completely remove/disassemble plunger. Place assembled cartridge/mix tip (with attached intra-oral tip, if applicable) through syringe side opening, sliding mix tip through barrel end with cross notches. Push cartridge/mix tip assembly into place. After loading, firmly pull on the end of the mix tip to make sure the loaded cartridge/mix tip is fully seated in the syringe.
6. Away from the patient field, depress plunger until both members fully engage cartridge. Continue to depress plunger until material flows from both ports into the mix tip. Bleed and discard a small amount of mixed material, then proceed immediately to clinical application. Mix tip may be gently curved to allow access to the tooth preparation if desired. To dispense material, apply a slow, steady pressure to plunger (using palm or thumb, as preferred). Excessive force is not necessary. If resistance is encountered, or if excessive force is required, remove the syringe from the patient field and check for plunger obstructions. If cartridge/mix tip obstruction is suspected, remove cartridge/mix tip assembly and replace, following steps above.

7. Allow material to set completely before attempting disassembly. To disassemble, withdraw the plunger, straighten mix tip if bent, and tap the mix tip against the counter to dislodge the cartridge/mix tip assembly. Properly dispose used cartridge/mix tip assembly.
8. Prepare the digit™ targeted delivery system syringe for subsequent reuse by following instructions below for cleaning and disinfection/sterilization.

IMPRESSION TECHNIQUES

1. **Dual Phase, one-step technique:** with a Rigid Tray or with Double Arch Technique

Technique Tip: It is highly recommended to utilize an assistant to assure simultaneous syringing and tray loading.

1.1 Select and/or prepare a suitable metal tray, a firm disposable tray or a custom tray. Use rigid trays of sufficient size to provide at least 2-3 mm thickness of impression material. Caulk First Bite® Impression Trays are indicated for the impression of one or two teeth prepared for crowns, inlays or onlays while simultaneously recording the opposing dentition and bite registration of the two arches.

1.2 Brush a thin layer of tray adhesive onto tray following tray adhesive manufacturer's *Directions for Use*. Adhesive should be used even with perforated trays.

1.3 If there is any debris or blood, clean the field with air water spray.

1.4. Remove excess water with suction. **Do not desiccate the tooth.** Dry and isolate with cotton rolls.

For syringe material:

1.5 Backfill the impression syringe with Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material or attach an intraoral tip to the end of the mix tip for direct intraoral syringing from a cartridge or prepare digit™ targeted delivery system handpiece and unit dose cartridge assembly.

For tray material:

- 1.6a Mix and place or dispense selected tray material (assistant) directly into impression tray (*see complete Directions for Use for selected tray material*). Do not layer the material or place it into rows. This technique ensures the least amount of air incorporation. **Technique tip:** With Regular Set wash and tray material, loaded tray must be inserted within 1'10" (one minute ten seconds) from time of first syringing of wash material intraorally.
 - 1.6b Remove cotton rolls. *It is suggested to simultaneously syringe (dentist) completely around the margins of the preparation/abutment, making sure to keep intraoral tip within impression material.* Inject syringe material into any existing anatomy and continue syringing around the preparation(s)/abutment(s) until the preparation/abutment(s) is completely covered with syringe material Do not blast with air. Aquasil Ultra Impression Material's hydrophilic properties improve adaptation to the tooth enabling better "wetting".
 - 1.7 Seat loaded tray. *Tray must be seated within 1'10" (one minute ten seconds) from time of first introduction of wash material intraorally.* Retain impression in position until firmly set. (**NOTE:** The material sets faster in the mouth than on the bench.) Check impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time (Regular Set) is 5 minutes from the start of mix. At 72°F (22°C) Aquasil Regular Set Impression Materials have a minimum work time of 2 minutes 15 seconds. At 98°F (37°C) Aquasil Regular Set Impression Materials have a maximum intraoral work time of 1'10". Variations in temperature from that stated above will affect work and set times. **Higher temperatures reduce work times and lower temperatures increase them.**
 - 1.8 Remove impression by pulling slowly to break seal, then snap out along the long axis of the tooth. Rinse impression under cold water and blow dry.
2. **Dual phase, Two-Step Technique**
Aquasil Putty, Reprosil® Putty and Reprosil® Quixx Putty are very high viscosity materials consisting of two pastes, which harden to form a base for final impressions. **Technique Tip:** It is highly recommended to utilize an assistant with a second dispensing gun to ensure simultaneous syringing and tray loading.

- 2.1 Prepare suitable tray, and complete the primary Aquasil Putty, Reprosil® Putty or Reprosil® Quixx Putty impression. See complete *Directions for Use* for Aquasil Putty, Reprosil® Putty or Reprosil® Quixx Putty. **NOTE:** Do not apply tray adhesive to the putty material when performing a two step technique.
- 2.2 If there is any debris or blood clean the field with air water spray.
- 2.3 Remove excess water with suction. **Do not desiccate the tooth.** Dry and isolate with cotton rolls.
- 2.4 Backfill the impression syringe with Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material or attach an intraoral tip to the end of the mix tip for direct intraoral syringing from a cartridge or prepare digit™ targeted delivery system handpiece and unit dose cartridge assembly.
- 2.5a Dispense Aquasil XLV or Aquasil Ultra LV Regular Set (assistant) directly into primary putty impression. Do not layer the material or place it into rows. This technique ensures the least amount of air incorporation.
- 2.5b Remove cotton rolls. *It is suggested to simultaneously syringe (dentist) completely around the margins of the preparation/abutment, making sure to keep intraoral tip within impression material.* Inject syringe material into any existing anatomy and continue syringing around the preparation(s)/abutment until the preparation/abutment is completely covered with syringe material.
- 2.6 Seat loaded tray. *Tray must be seated within 1'10" (one minute ten seconds) from time of introduction of wash material intraorally.* Retain impression in position until firmly set. (**NOTE:** The material sets faster in the mouth than on the bench.) Check impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time (Regular Set) is 5 minutes from the start of mix. At 72°F (22°C) Aquasil Regular Set Impression Materials have a minimum work time of 2 minutes 15 seconds. At 98°F (37°C) Aquasil Regular Set Impression Materials have a maximum intraoral work time of 1'10". Variations in temperature from that stated above will affect work and set times. **Higher temperatures reduce work times and lower temperatures increase them.**
- 2.7 Remove impression by pulling slowly to break seal, then snap out along the long axis of the tooth. Rinse impression

under cold water and air dry.

3. **Edentulous/Overdenture Impression Technique**

- 3.1 Prepare a suitable tray with proper border extensions. Higher viscosities of Aquasil Ultra or Reprosil® material may be used for border molding. See complete *Directions for Use* for selected border molding material.
- 3.2 Brush a thin layer of tray adhesive onto tray extending to the borders of the tray, following tray adhesive manufacturer's *Directions for Use*.
- 3.3 The final impression is made with Aquasil Ultra XLV or Aquasil Ultra LV Impression Material. Dispense the Aquasil Ultra XLV or Aquasil Ultra LV material onto impression tray. Do not layer the material or place it into rows. This technique ensures the least amount of air incorporation.
- 3.4 Seat loaded tray. Retain impression in position until firmly set. **(NOTE:** The material sets faster in the mouth than on the bench.) Check impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time (Regular Set) is 5 minutes from the start of mix. At 72°F (22°C) Aquasil Ultra Regular Set Impression Materials have a minimum work time of 2 minutes 15 seconds. Variations in temperature from that stated above will affect work and set times. **Higher temperatures reduce work times and lower temperatures increase them.**
- 3.5 Remove impression by pulling slowly to break seal, then snap out of the mouth. Rinse impression under cold water and air dry.

4. **Duplication of Models**

- 4.1 Select and/or prepare a suitable metal tray, a firm disposable tray or a custom tray. Use rigid trays of sufficient size to provide at least 2-3 mm thickness of impression material.
- 4.2 Brush a thin layer of tray adhesive onto tray following tray adhesive manufacturer's *Directions for Use*.
- 4.3 Load the impression tray with Aquasil Ultra Monophase, Heavy or Rigid Regular Set tray material. Do not layer the

material or place it into rows. This technique ensures the least amount of air incorporation.

- 4.4 It is suggested to simultaneously syringe, using Aquasil Ultra XLV or Aquasil Ultra LV Regular Set Impression Material, completely around the model teeth, making sure to keep tip within impression material. Inject syringe material into any existing anatomy and continue syringing around the teeth until they are completely covered with syringe material.
- 4.5 Seat model into the loaded tray. Retain impression in position until firmly set. **(NOTE:** The material sets slower on the bench than in the mouth.) Check impression to be sure it is firm, resilient, and non-tacky before removal. The Minimum Removal Time (Regular Set) is 5 minutes from the start of mix. At 72°F (22°C) Aquasil Ultra Regular Set Impression Materials have a minimum work time of 2 minutes and 15 seconds. Variations in temperature from that stated above will affect work and set times. **Higher temperatures reduce work times and lower temperatures increase them.**
- 4.6 Remove model from impression by pulling slowly to break seal.

CLEANING

DISPENSER AND CARTRIDGE CLEANING INSTRUCTIONS

The cartridge dispensing guns and digit™ targeted delivery system syringe may be cleaned by scrubbing with hot water and soap or detergent. The digit™ targeted delivery system syringe plunger should be fully retracted. Excess material may be cleaned with alcohol-moistened gauze. Do not attempt to clean, disinfect or re-use unit dose cartridge or tips. Properly dispose used tips and unit dose cartridges.

DISINFECTION

DISPENSER AND CARTRIDGE DISINFECTION INSTRUCTIONS

Dispensers, cartridges or used tips, if left in place on cartridges for storage, exposed to splatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant.

Acceptable disinfectants are those that are EPA-registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Water-based disinfectant solutions are preferred. Use impregnated wipes not sprays. **NOTE:** Wipe cartridge GENTLY, vigorous wiping may destroy the label.

Some phenolic based agents and iodophor based products may cause surface staining. Agents containing organic solvents, such as alcohol may tend to dissolve the plastic dispenser and plunger. The disinfectant manufacturer's directions should be followed properly for optimum results. **Technique Tip:** Unset Aquasil Impression Material can be removed from clothing with dry cleaning solvent.

DISINFECTION AND/OR STERILIZATION

DIGIT™ TARGETED DELIVERY SYSTEM SYRINGE

Following cleaning as outlined above, digit™ targeted delivery system syringe exposed to splatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be steam autoclaved following autoclave manufacturer's recommendations.

Alternatively, the digit™ targeted delivery system syringe may be disinfected with a hospital-level disinfectant as outlined above. Water-based disinfectant solutions are preferred. Some phenolic-based agents and iodophor-based products may cause surface staining. Agents containing organic solvents, such as alcohol may tend to dissolve the plastic barrel and plunger. The disinfectant manufacturer's directions should be followed properly for optimum results. **NOTE:** As with any plastic instrument, the digit™ targeted delivery system syringe may weaken over time.

To reassemble, insert plunger into syringe barrel, and press components together. Prior to each use check to make sure that the digit™ targeted delivery system syringe plunger is fully engaged and in good working order.

DISINFECTION OF THE IMPRESSION MATERIAL

The impression should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are those that are EPA-registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide, and dual or synergized quaternary ammoniums are approved disinfectants. To disinfect polyvinyl siloxane material, thoroughly soak by spraying or immerse the impression in any recommended hospital level disinfectant except neutral glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Water-based disinfectant solutions are preferred. To avoid surface porosity in stone, delay pour at least 15 minutes after disinfection; for epoxy dies delay pour 60 minutes.

CASTING AND ELECTROPLATING

The impression should be removed from the disinfectant, rinsed with water and exposed to air to dry for at least the same amount of time it was exposed to the disinfectant solution before it is poured. The impression may be poured up to fourteen days. Do not store impression in direct sunlight. If the impression is to be shipped, use suitable packaging to preclude distortion. The material is compatible with a range of die stones. Impressions may be silver or copper-plated.

LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: "YYYY/MM."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number on the cartridge
 - Expiration date

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©2005 DENTSPLY International Inc. All Rights Reserved. Printed in U.S.A.
Manufactured in the U.S.A. by:
DENTSPLY DeTrey GmbH