CLEARFIL MAJESTY™ Esthetic

CE
0197 93/42/EEC

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL MAJESTY Esthetic is a light-cure, radiopaque restorative composite resin which provides accurate color matching, high polish ability and excellent physical properties, making it ideal for both anterior and posterior restorations. It is formulated with optimal viscosity assuring easy handling and placement.

II. INDICATIONS

CLEARFIL MAJESTY Esthetic is indicated for the following restorative applications:

- Direct restorations for anterior and posterior teeth (Class I – V cavities)
- Direct veneers
- Correction of tooth position and tooth shape (e.g. diastema closure, dwarfed tooth, etc.)
- Intraoral repairs of fractured crowns / bridges

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers.

IV. INCOMPATIBILITIES

Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.

V. PRECAUTIONS

1. *Safety precautions*

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomer.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
 - Immediately rinse the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin and the oral mucosa>
 - Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and rinse with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the dental curing light when curing the product.
- To prevent cross infection, avoid sharing the same paste among different patients. Dispense the paste on a mixing paper before use. Do not try to apply the paste directly from the syringe onto the patient's tooth.

2. *Handling and manipulation precautions*

- Do not use the product in conjunction with the other composite resin. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- If the adherent surface extends to uncut enamel, apply a phosphoric acid etching agent (e.g. K-ETCHANT GEL) and leave it in place for 10 seconds before washing and drying.
- Amalgam or other lining material remaining in the cavity will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the cavity.
- When light curing the product, note the light curing depth in this Instructions for Use.
- The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light-evaluating device at appropriate intervals.
- The paste should be used as soon as practicable after being dispensed from the syringe. If the paste is to be left for a while before use, it should be covered with a light blocking plate.

- The product should be returned to room temperature before dispensing. If it has been taken out of a refrigerator. Failure to do this might cause breakage of the syringe.
- The use of the product is restricted to licensed dentist.

3. *Storage precautions*

- Do not use the product after the expiration date. Note the expiration date on the outside of package.
- The product should be stored at 2 - 25°C / 36 - 77°F when not in use.
- Keep away from extreme heat, direct sunlight or flame.
- After the paste has been dispensed, the syringe should be capped securely as soon as possible to prevent the paste inside the syringe from being hardened by ambient light, and to keep foreign matter from entering the syringe.
- The product must be stored in proper places, where only dental practitioners are able to handle.

VI. SHADES SYSTEM AND COMPONENTS

1. *Components*

- CLEARFIL MAJESTY Esthetic is available in 17 shades;
 - A1, A2, A3, A3.5, A4, B2, B3, C3, HO, E, XL, OC, T, AM, OA2, OA3, OA4
- Shade guide

2. *Ingredients*

- Principal ingredients:
 - Silanated barium glass powder
 - Pre-polymerized organic filler
 - Bis-phenol A diglycidylmethacrylate (Bis-GMA)
 - Hydrophobic aromatic dimethacrylate
 - dl-Camphorquinone

VII. CLINICAL PROCEDURES

Method. A

- Direct restorations for anterior and posterior teeth (Class I – V cavities)
- Direct veneers
- Correction of tooth position and tooth shape (e.g. diastema closure, dwarfed tooth, etc.)

1. Cleaning tooth structure

Be sure the cavity is adequately cleaned, An adequately cleaned cavity assures maximum adhesive performance.

2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

3. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

5. Tooth surface treatment and bonding

Tooth surface treatment and bonding according to Instructions for Use of the Bonding system used (e.g. CLEARFIL TRI-S BOND).

6. Placement and light curing of CLEARFIL MAJESTY Esthetic

6-1. If necessary, the application of a flowable composite resin (e.g. CLEARFIL LINER BOND 2 PROTECT LINER F) to the cavity as an initial layer is optional (Refer to its' Instructions for Use). This layer should be polymerized prior to placing the first layer of CLEARFIL MAJESTY Esthetic.

6-2. Select an appropriate shade of the product, using Shade guide, corresponding to the shade of the tooth you are going to treat. You may combine two or more shades to achieve a better match to the shade of the tooth to be treated.

6-3. Place the chosen shade product into the cavity, and light-cure the resin with a dental curing unit* following the table.

Table: Relation between curing time and depth of cure for each dental curing unit.

Type	Curing time (Sec)	Depth of cure (mm)			
		T, AM	A1,A2,A3,A3.5,A4,B2,B3,C3,HO,E,XL,OC	OA2, OA3, OA4	
Conventional halogen*	20	2.0	1.5	1.0	
Fast halogen*	10				
Plasma arc*	5				
LED*	20				

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity ²⁾ of 300 - 550 mW/cm ² in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity ²⁾ of more than 550 mW/cm ² in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity ²⁾ of more than 2000 mW/cm ² in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm ² in wavelength range from 400 - 430 nm.
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

1) Peak of emission spectrum: 450 - 480 nm.

2) Wavelength distribution and light intensity values measured with a spectroradiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp.

[NOTE]

- If two or more shades are selected, the paste of one shade should be applied and then light cured before the paste of another shade is applied and light cured. After all the pastes have been layered, the entire bulk of the product should be completely light cured again.
- If necessary, a thin layer of opaque resin (e.g. CLEARFIL ST OPAQUER) can be applied and light-cured according to its' Instructions for Use.

7. Finishing

Contour the restoration and adjust the occlusion using a fine diamond point. Polish with silicon rubber points or polishing discs.

Method. B

- Intraoral repairs of fractured crowns / bridges

- Preparation of fractured surfaces
 - a. Facing material surface
 - If necessary, remove a layer of the fractured surfaces and place a bevel at the marginal area using suitable instruments (e.g. diamond point).
 - b. Metal surface
 - If the metal surface is not surface treated, roughen it by using sandblasting with alumina powder or using a diamond point. Rinse the surface with water and then dry.
- Acid etching of facing material
 - Apply a phosphoric acid etching agent (e.g. K-ETCHANT GEL) to the facing material and leave it in place for 5 seconds before washing and drying.
- Precious metal surface
 - If precious metal is used, apply a metal-adhesive primer (e.g. ALLOY PRIMER) according to the Instructions for Use of the product used.
- Silane treatment / bonding
 - 4-1. Apply a silane-coupling agent (e.g. the mixture of CLEAFIL TRI-S BOND and CLEARFIL PORCELAIN BOND ACTIVATOR) to the facing material according to the Instructions for Use of the product used.
 - 4-2. If the adherent surface extends to the tooth, apply a bonding agent (e.g. CLEARFIL SE BOND, CLEARFIL TRI-S BOND) according to the Instructions for Use of the product used.

- Placement and light curing of CLEARFIL MAJESTY Esthetic
 - 5-1. If necessary, the application of a flowable composite resin (e.g. CLEARFIL LINER BOND 2 PROTECT LINER F) to the cavity as an initial layer is optional (Refer to its' Instructions for Use). This layer should be polymerized prior to placing the first layer of CLEARFIL MAJESTY Esthetic.

5-2. Using Shade guide, select an appropriate shade of the product corresponding to the fractured crown / bridge to be repaired. Two or more shades can be selected and combined in order to give a better match to the shade of the fractured crown / bridge.

5-3. Load CLEARFIL MAJESTY Esthetic of the selected shade to the cavity and light-cure it with a dental curing unit for curing time specified in table in **Method A** 6-3. Considering the depth of cure, apply the resin in layers and polymerize each layer, if

necessary. If two or more shades are to be applied, apply them in layers respectively, either polymerizing all the layers collectively, or each layer individually, depending on the total depth of cure.

[NOTE]

Use opaque resin (e.g. CLEARFIL ST OPAQUER) for metal surface to prevent metal shine through.

6. Finishing

Contour the restoration and adjust the occlusion using a fine diamond point. Polish with silicon rubber points or polishing discs.

[WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL MAJESTY is a trademark of KURARAY CO., LTD.
CLEARFIL MAJESTY ist ein eingetragenes Warenzeichen der KURARAY CO., LTD.
KURARAY EUROPE GmbH