DUAL CURE DENTAL ADHESIVE SYSTEM **PANAVIA™ F 2.0**

CE
0197
93/42/EEC

GB INSTRUCTIONS FOR USE

I. INTRODUCTION

PANAVIA F 2.0 is a dual-curing (chemical and light) resin based cement system for metal, composite, and silanated porcelain restorations. PANAVIA F 2.0 consists of ED PRIMER II, PANAVIA F 2.0 paste and OXYGUARD II. ED PRIMER II contains HEMA and 5-NMSA as well as MDP and consists of liquid A and liquid B. PANAVIA F 2.0 Paste releases fluoride. ALLOY PRIMER improves the bond strength for precious alloy and PANAVIA F 2.0 Paste. Tin plating is not required.

II. INDICATIONS

PANAVIA F 2.0 is indicated for the following applications:

- [1] Cementation of metal crowns and bridges, inlays and onlays.
- [2] Cementation of porcelain crowns, inlays, onlays and veneers.
- [3] Cementation of composite resin crowns, inlays, and onlays.
- [4] Cementation of adhesive bridges.
- [5] Cementation of endodontic cores and prefabricated posts.
- [6] Amalgam bonding.

[NOTE]

Use cement shades appropriate to individual cases. The shades of adhesive cement and applicable cases

Cement shade	TC, Light	White	Opaque
Metal inlays and onlays; metal crowns and bridges	○	○	○
Porcelain or composite inlays, onlays and veneers;	○	△	△
Preformed posts and cast metal cores	○	○	○
Adhesion bridges and splints anterior	△	○	○
posterior	○	○	○
Bonded amalgam restorations	○	○	○

○=RECOMMENDED ○=ADEQUATE △=NOT RECOMMENDED

III. CONTRAINDICATIONS

[1] Patients with allergy to methacrylate monomers.
[2] Patients with allergy to acetone (ALLOY PRIMER)

IV. SIDE EFFECTS

The mucous membrane may turn whitish when contacted by primer due to protein coagulation. This is a temporary phenomenon that will disappear in a few days.

V. INCOMPATIBILITY

[1] Do not use eugenol containing materials for pulp protection and temporary sealing since they retard the curing process.
[2] Do not use hemostatics containing ferric compounds, since they may impair adhesion and the remaining ferric ion may cause discoloration in the margin or surrounding gingiva.
[3] Do not use PANAVIA F 2.0 paste with PANAVIA F paste, or the mixed paste will not cure completely by light curing.

VI. PRECAUTIONS

1. Safety precautions

If any hypersensitivity such as a rash or dermatitis occurs due to contact with the product, discontinue the use of the product and consult a physician.
2. In order to prevent the occurrence of hypersensitivity due to contact with methacrylate monomers, acetone or acid, wear gloves or take other appropriate protective measures.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using product, cover the patients eyes with a towel to protect the patient eyes from splashing material.
4. If the product comes in contact with the human body, take the following actions:
-<i>If the product gets in the eye:</i> Immediately wash the eye with copious amounts of water and consult a physician.
-<i>If the product comes in contact with the skin:</i> immediately wipe it off with a cotton pledge moistened with alcohol or gauze and wash with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. If a patient or operator feels sick from inhaling the acetone contained in the product, allow them to rest and breathe fresh air.
7. Any actual or near pulp exposure area should be covered with a hard setting calcium hydroxide material. Do not use eugenol materials for pulp protection.
8. When using with preformed stainless posts, the post should not contact metal restorations. Cover the post with composite resin.
9. Avoid looking directly at the curing light when curing the product.
10. Keep the product out of reach of children.

2. Handling and manipulation precautions

1. PANAVIA F 2.0 polymerizes by a dental curing light (irradiate wave length : 400-515nm). Use the light blocking plate to avoid exposing the material to an operating light or natural light (sunlight from windows).
2. Make sure the disposable nozzle or disposable brush tip is securely attached to prevent the patient from swallowing them.
3. After placing the restoration on the tooth, the cement could cure due to the operatory light. Use caution not to let the operatory light come too close to the patient.

3. Storage precautions

1. Do not use after the expiration date. Note expiration date on the outside of package.
2. ALLOY PRIMER is flammable. Keep away from flame.
3. The product should be refrigerated (2-8°C/36-46°F) when not in use; bring it to room temperature before using.
4. ALLOY PRIMER should be stored at 2-25°C/36-77°F when not in use.
5. Keep away from extreme heat or direct sunlight.
6. The bottle or syringe cap should be replaced as soon as the resin has been dispensed from the bottle or syringe. This prevents evaporation of volatile ingredients.

VII. KIT COMPONENTS

PANAVIA F 2.0 consists of the following components. They are available as a set or single components.

- 1) ED PRIMER II (adhesive primer) Liquid A 4ml, Liquid B 4ml
* There is no need for washing after application and the primer should be left in place.
- 2) PANAVIA F 2.0 paste A paste 2.3ml, B paste 2.3ml
* B paste is available in four shades: TC, Light, White and Opaque.
- 3) ALLOY PRIMER 1 ml (Refill bottle holds 5ml)
* This metal adhesive primer provides high bond strength to both precious and base metals.
* It enhances the bond strength of PANAVIA F 2.0 to precious metals.
- 4) OXYGUARD II 6 ml
* This Oxygen-blocking agent allows polymerization of the PANAVIA F 2.0 Paste.
- 5) Accessories Mixing plate, spatula, mixing paper, small brush holder, disposable brush tip, disposable nozzle, light-blocking plate

VIII. RELATED PRODUCTS

The following products are necessary for specific procedures.

- 1) CLEARFIL PORCELAIN BOND ACTIVATOR 4 ml
* This product contains a silane coupling agent. Mixing it with CLEARFIL SE BOND PRIMER or CLEARFIL LINER BOND 2V PRIMER or CLEARFIL NEW BOND or CLEARFIL PHOTO BOND improves the bond strength to porcelain or cured composite.
- 2) K ETCHANT GEL 6 ml
* This phosphoric acid solution is used for pretreatment of uncut enamel and porcelain.

IX. CLINICAL PROCEDURES

1. Basic procedure (use of adhesive cement)

[NOTE]
Use the mixed paste as soon as possible after dispensing and mixing.

- [1] Dispensing the pastes**
1. Align marking on the nut with the reference line on the plunger and turn the syringe to dispense the necessary amount of A paste minimum turn of the syringe should be half a turn.
2. Equal amounts of A paste and B paste should be dispensed.
3. The amount of paste dispensed at last rotation of the syringe could be inaccurate. Therefore, discard the syringe before using the last portion.
4. The necessary amount of paste for a typical applications is:

Number of rotations of syringe	Applications
Half a turn	Inlays and onlays
1 turn	Crowns

[NOTE]
1. If the paste is dispensed by turning the plunger a quarter of a turn, the performance of the product could be impaired when the paste hardens.
2. If not used immediately, it should be covered with a light-blocking plate.

[2] Mixing A paste and B paste
Mix sufficient A paste and B paste on the mixing plate for 20 seconds. Be sure there is no water mist on the mixing plate or spatula before using them; the presence of water could shorten the working time of the mixed paste.

[3] Preparing PANAVIA F 2.0 paste
Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[4] Cementation

1. Applying mixed paste to the restoration

Apply the mixed paste to the restoration.

[5] Cleaning of tooth structure

▼ Clean and prepare the cavity in the normal manner for an amalgam restoration

[6] Treatment of tooth surface

▼ Treat the entire cavity surface with ED PRIMER II

[7] Preparing of adhesive cement

▼ Dispense the necessary amount of each paste according to the case and mix the two paste

[8] Placing the amalgam

▼ Apply the mixed paste to the entire cavity surface and place the triturated amalgam into the abutment or cavity

[9] Removing excess cement

▼ Remove excess cement with an explorer or small scaler. The restoration can then be finished and polished with pumice and water.

[10] Curing cement margin

Cure the mixed paste along the cement margin, using either of the following two methods.

① Light curing method

When it is possible to light cure adhesive cement along the cement margin, such as inlays and onlays, light cure each section of the cement margin for 20 seconds by conventional halogen curing lights (>250mW/cm²) or LED curing lights (>160mW/cm²). If plasma arc curing lights (>2000mW/cm²) or fast halogen curing lights (>550mW/cm²) are used, each section of the cement margin can be cured for 5 seconds.

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[CAUTION]

If the adherend surface is contaminated with saliva or blood after ALLOY PRIMER is applied, clean the adherend surface with a cotton pledge moistened with alcohol and apply ALLOY PRIMER again.

[3] Applying ED PRIMER II

Apply ED PRIMER II to the entire tooth surface (enamel and dentin) of the adherend surface or metal or composite resin abutment tooth with a disposable brush tip or sponge and leave it in place for 30 seconds.

[4] Drying

Using a sponge or paper point, remove excess primer to prevent the formation of a pool of the primer in the root canal or at the corners of the cavity. Dry the primer completely with gentle air flow. Remember that a pool of the primer will cause quick polymerization of the adhesive cement. It is necessary to fit the post as quickly as possible.

[5] Spreading excess cement

Using a small brush, spread excess cement over the remaining crown and post head.

[6] Curing adhesive cement

Light cure the adhesive cement on the remaining crown and post head for 20 seconds by conventional halogen curing lights (>250mW/cm²) or LED curing lights (>160mW/cm²). If plasma arc curing lights (>2000mW/cm²) or fast halogen curing lights (>550mW/cm²) are used, each cement margin has to be cured for 5 seconds.

[7] Finishing

Never use a lentalo spiral for loading the adhesive cement into the root canal. If the adhesive cement is loaded into the root canal using a composite resin syringe, the polymerization of the cement is accelerated. It is necessary to fit the post as quickly as possible.

[8] Sealing the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrainment of air bubbles.

[9] Adhering surface treatment

After treating the restoration surface, proceed to cementation quickly.

[10] Adherend surface treatment

1) Cleaning cavity or abutment tooth surface

2) Applying temporary sealing agent or temporary cementation material from the adherend surface.

[11] Preparing PANAVIA F 2.0 paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[12] Sealing the post

1. Apply the mixed paste to the post.

2. Seal the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrainment of air bubbles.

[13] Adhering surface treatment

After treating the restoration surface, proceed to cementation quickly.

[14] Preparing PANAVIA F 2.0 paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[15] Sealing the post

1. Apply the mixed paste to the post.

2. Seal the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrainment of air bubbles.

[16] Adhering surface treatment

After treating the restoration surface, proceed to cementation quickly.

[17] Preparing PANAVIA F 2.0 paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[18] Sealing the post

1. Apply the mixed paste to the post.

2. Seal the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrainment of air bubbles.

[19] Adhering surface treatment

After treating the restoration surface, proceed to cementation quickly.

[20] Preparing PANAVIA F 2.0 paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[21] Sealing the post

1. Apply the mixed paste to the post.

2. Seal the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrainment of air bubbles.

[22] Adhering surface treatment

After treating the restoration surface, proceed to cementation quickly.

[23] Preparing PANAVIA F 2.0 paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

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