

Filtek™ Supreme XT Flowable Restorative

A new benchmark for flowable composites delivering high strength *and* low wear.



- 3M ESPE's unique nanotechnology results in better esthetics
- Better gloss retention than leading flowables
- Flow-on-demand handling that dentists prefer
- Better in vitro wear than leading flowables and comparable to some traditional microhybrids
- Easy to use
- Excellent radiopacity
- Available in 12 shades including Opaque A3 which makes a great liner/base

More dentists are discovering the advantages of using a “flowable composite” in restorative procedures.

Flowable composites contain less filler creating a lower viscosity – making them “more fluid” than traditional composites. That means they can flow more easily and cover areas more effectively than other composites.


3M ESPE applied our in-depth knowledge of nanotechnology to create **Filtek Supreme XT Flowable Restorative** – with many of the desirable attributes of Filtek™ Supreme XT Universal Restorative.

- **Flow-on-demand handling dentists prefer:** Like Filtek™ Flow Flowable Restorative, it flows under pressure yet holds its shape and stays in place prior to light curing. No oozing, slumping or running.

- **Excellent esthetic properties** including polish, polish retention and color match, for a higher level of patient satisfaction.

- **Available in 12 shades** that match Filtek™ Supreme XT Universal Restorative.

It looks great, and wears well, too! *See back for details.*




Test data proves Filtek Supreme XT Flowable Restorative retains its polish better than other leading flowables. SOURCE: 3M ESPE Internal Data

Filtek™ Supreme XT Flowable Restorative wears better than other leading flowables.

Strength derived from 3M ESPE's nanotechnology gives you real confidence using a flowable in various restorative procedures.

- **Excellent compressive and tensile strength and optimal wear resistance** make it ideal for a variety of indications, including as a base/liner, pit & fissure sealant, and Class III and V restorations.
- **Ideally suited for Minimally Invasive Dentistry (MID)**



In vitro wear of Filtek Supreme XT Flowable Restorative is better than leading flowables and comparable to some traditional microhybrids.
SOURCE: 3M ESPE Internal Data

Ordering Information

Filtek Supreme XT Flowable Restorative

Item Number	Product Information
3913A1	2 Syringes, 2 g each, Shade A1, 20 Single-Use Application Tips
3913A2	2 Syringes, 2 g each, Shade A2, 20 Single-Use Application Tips
3913A3	2 Syringes, 2 g each, Shade A3, 20 Single-Use Application Tips
3913A3.5	2 Syringes, 2 g each, Shade A3.5, 20 Single-Use Application Tips
3913A4	2 Syringes, 2 g each, Shade A4, 20 Single-Use Application Tips
3913B1	2 Syringes, 2 g each, Shade B1, 20 Single-Use Application Tips
3913B2	2 Syringes, 2 g each, Shade B2, 20 Single-Use Application Tips
3913C2	2 Syringes, 2 g each, Shade C2, 20 Single-Use Application Tips
3913D2	2 Syringes, 2 g each, Shade D2, 20 Single-Use Application Tips
3913W	2 Syringes, 2 g each, Shade W, 20 Single-Use Application Tips
3913XW	2 Syringes, 2 g each, Shade XW, 20 Single-Use Application Tips
3913OA3	2 Syringes, 2 g each, Shade OA3, 20 Single-Use Application Tips
3700T	20 Single-Use Application Tips

Introductory Kit

3913	5 Syringes, 2 g each, Shades A2, A3, A3.5, W, OA3, 60 Single-Use Application Tips
------	---



3M ESPE AG · ESPE Platz
82229 Seefeld · Germany
E-Mail: info3mespe@mmm.com
Internet: http://www.3mespe.com

3M, ESPE and Filtek are trademarks of 3M or 3M ESPE AG. Tetric, TPH3, EsthetX, Flow-It, Revolution and Point 4 are not trademarks of 3M or 3M ESPE AG.

© 3M 2005. All rights reserved.