

Ketac™ Silver Aplicap™/Maxicap™ Glass Ionomer Filling Material



Ketac Silver restorative is available in a convenient Aplicap or Maxicap delivery for core build-ups, small posterior and primary tooth restorations.

When patient satisfaction is your #1 priority...use the #1 metal-reinforced glass ionomer solution.

Ketac Silver glass ionomer restorative is a proven alternative to amalgam for core build-ups.

A patented 3M ESPE process sinters silver particles into the glass for excellent thermal and electrical insulation – thus reducing post-operative sensitivity.

Moisture tolerance, long-term fluoride release, and adhesion to enamel and dentin give you added confidence of sealed, lasting restorations.



Ketac Silver material is a fluoride-releasing alternative to amalgam for large fillings.

Clinical advantages:

- Ketac Silver is biocompatible. This makes it particularly suitable for core build-ups.
- An excellent insulator, which protects the pulp and helps eliminate patient sensitivity.
- No need for undercutting – so you conserve tooth structure.
- An easy-to-mix hand mix formula and two capsule delivery systems makes it versatile in any application.
- Less shrinkage than traditional composites – leads to better marginal integrity.

Indications:

- Core build-ups
- Primary teeth fillings
- Small fillings in posterior teeth
- Temporaries
- Base under amalgams













Ketac™

- Contains no amalgam.
- High sustained fluoride release.
- Chemical bond to enamel and dentin provides excellent marginal integrity.
- Radiopaque.
- Aplicap and Maxicap delivery systems provide for consistent mixing and easy dispensing.
- Maxicap provides four times the volume of Aplicap. Designed for larger scale restorations and core build-ups.

Glass Ionomer solutions from 3M ESPE produce the results you depend on.

3M ESPE is a trusted innovator of glass ionomers, with clinically proven solutions for your procedural needs. Use this quick reference – based on more than 30 years of experience – to select the material best suited for your patient.

Indications										
	Small Class I	Class III	Class V	Minimal Invasive Dentistry (MID)	Primary Teeth	Base Under Composite	Temporary Treatment	Core Build-up	Aplicap™ Delivery	Maxicap™ Delivery
CONVENTIONAL										
Ketac Molar Quick	●		■	●	●	●	●	●	★	
Ketac Molar	●		■	●	●	●	●	●	★	
Ketac Fil Plus	●	●	●	●	●			●	★	
RESIN MODIFIED										
Photac Fil	●	●	●	●	●		●		★	
Vitremer™	●	●	●	●	●	●	●	●		
METAL REINFORCED										
Ketac Silver	●		■		●		●	●	★	★

■ = non esthetic

Ordering Information



Ketac Silver Maxicap Introductory Kit (56204)

3M™ ESPE™ RotoMix™ Capsule Mixing Unit (76310)



Item Number	Product Information
37801	Ketac Silver Handmix Introductory Kit Kit contains: 1 – 25g Powder; 1 – 12ml Liquid; Mixing Pad
37811	Ketac Silver Handmix Refills 1 – 25g Powder
37821	Ketac Silver Handmix Refills 1 – 12ml Liquid
37007	Ketac Silver Aplicap Introductory Kit Kit contains: 50 Aplicap Capsules; 1 Activator; 1 Applier; Ketac Glaze Light-Cured Varnish; Ketac Conditioner Dentin Pretreatment Solution; Mixing Pad
37011	Ketac Silver Aplicap Refills 50 Capsules
56204	Ketac Silver Maxicap Introductory Kit Kit contains: 50 Maxicap Capsules; 1 Activator; 1 Applier; Ketac Glaze Light-Cured Varnish; Ketac Conditioner Dentin Pretreatment Solution; Mixing Pad
56211	Ketac Silver Maxicap Refills 50 Capsules
37461	Ketac™ Conditioner 1 – 10ml Liquid Pack, 25 Tips
37191	Ketac™ Glaze 1 – 2.5g Bottle
76310	RotoMix Capsule Mixing Unit (120V/60Hz) Simple touch-pad, programmable operation. Delivers a completely homogenous mix. Designed for mixing 3M ESPE Aplicap and Maxicap capsules but is suitable for all standard dental capsule mixing materials.

For More Information

For more information on 3M ESPE Glass Ionomer Restorative Materials, visit our Web site or call our technical hotline. Products are available for purchase through authorized 3M ESPE distributors.

Technical Hotline: 1-800-634-2249 www.3MESPE.com




Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
USA

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-265-1840 ext. 6229



Printed in USA
© 3M 2004. All rights reserved.

3M, ESPE, Ketac, Photac, Vitremer, Aplicap, Maxicap and RotoMix are trademarks of 3M or 3M ESPE AG.
70-2009-3715-2