

Ketac™ Fil Plus Aplicap™

Glass Ionomer Filling Material



Ketac Fil Plus Aplicap capsule delivery provides a time-saving, precise mixing ratio for outstanding, reproducible results.

The esthetic glass ionomer solution you can count on for restoring a confident smile.

Choose Ketac Fil Plus material to closely match natural tooth structure.


Ketac Fil Plus filling material is a conventional glass ionomer that helps you achieve virtually invisible anterior restorations.

There are eight shades to choose from, making it easy to match natural tooth color.

And, with the chemical bond to enamel and dentin, less tooth structure needs to be removed during preparation – which supports a healthy smile.

Ketac Fil Plus mirrors natural tooth structure.

Coefficient of Thermal Expansion



Source: Internal data 3M ESPE AG

Clinical advantages:

- Ideal for pediatric and geriatric patients.
- Can be used to treat lost or cervical restorations.
- Less shrinkage than traditional composites – leads to better marginal integrity.
- Unique capsule design assures proper activation – so time-wasting, costly “duds” are not a problem.

Indications:

- Class III and V restorations
- Wedge-shaped defects
- Small Class I fillings
- Core build-ups
- Minimally Invasive Dentistry (MID)











Ketac™



- Chemical bond to enamel and dentin guarantees a tight marginal seal.
- Tooth-like coefficient of thermal expansion.
- Aplicap capsule delivery for consistent handling characteristics and reproducible results.
- Minimal preparation.
- High fluoride release.
- Radiopaque.
- Available in eight esthetic shades.

Glass Ionomer solutions from 3M ESPE produce the results you depend on.

3M ESPE is a trusted innovator of glass ionomers, with clinically proven solutions for your procedural needs. Use this quick reference – based on more than 30 years of experience – to select the material best suited for your patient.

Indications										
	Small Class I	Class III	Class V	Minimal Invasive Dentistry (MID)	Primary Teeth	Base Under Composite	Temporary Treatment	Core Build-up	Aplicap™ Delivery	Maxicap™ Delivery
CONVENTIONAL										
Ketac Molar Quick	●		■	●	●	●	●	●	★	
Ketac Molar	●		■		●	●	●	●	★	
Ketac Fil Plus	●	●	●	●	●	●	●	●	★	
RESIN MODIFIED										
Photac Fil	●	●	●	●	●		●		★	
Vitremer™	●	●	●	●	●	●	●	●		
METAL REINFORCED										
Ketac Silver	●		■		●		●	●	★	★

■ = non esthetic

Ordering Information



Ketac™ Fil Plus Aplicap™ Introductory Kit (55001)

3M™ ESPE™ RotoMix™ Capsule Mixing Unit (76310)



Item Number	Product Information
55001	Ketac Fil Plus Aplicap Introductory Kit Kit contains: Activator; Applier; 50 Assorted Aplicap Capsules (10 each of A1, A2, A3, A3.5 and A4)
55011	Ketac Fil Plus Aplicap 50 – Capsules, Assorted
55021	Ketac Fil Plus Aplicap 50 – A1 Shade Capsules
55031	Ketac Fil Plus Aplicap 50 – A2 Shade Capsules
55041	Ketac Fil Plus Aplicap 50 – A3 Shade Capsules
55051	Ketac Fil Plus Aplicap 50 – A3.5 Shade Capsules
55061	Ketac Fil Plus Aplicap 50 – A4 Shade Capsules
55071	Ketac Fil Plus Aplicap 50 – C2 Shade Capsules
55081	Ketac Fil Plus Aplicap 50 – C3 Shade Capsules
55091	Ketac Fil Plus Aplicap 50 – DGO Shade Capsules
37461	Ketac™ Conditioner 1– 10ml Liquid Pack, 25 Tips
37191	Ketac™ Glaze 1– 2.5g Bottle
76310	RotoMix Capsule Mixing Unit (120V/60Hz) Simple touch-pad, programmable operation. Delivers a completely homogenous mix. Designed for mixing 3M ESPE Aplicap and Maxicap capsules but is suitable for all standard dental capsule mixing materials.

For More Information

For more information on 3M ESPE Glass Ionomer Restorative Materials, visit our Web site or call our technical hotline. Products are available for purchase through authorized 3M ESPE distributors.

Technical Hotline: 1-800-634-2249 www.3MESPE.com




Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
USA

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-265-1840 ext. 6229



Printed in USA
© 3M 2004. All rights reserved.

3M, ESPE, Ketac, Photac, Vitremer, Aplicap, Maxicap and RotoMix are trademarks of 3M or 3M ESPE AG.
70-2009-3714-5