



ISO 9001 Cert Num 4637/0 - 4638/0 - IT-53581 - IT53582



“ Technical Data”

MANUFACTURER:	BY DENTAL srl Pistoia – Italy www.bydental.com
MODEL:	Sky Light Export
DIMENSION:	APPROX h90
MATERIAL: Head Pole Base	Bay Bland Iron Iron
WEIGHT:	∞4 Kg
LED:	12 High Intensity
POWER OF THE LED:	8 W
WAVELENGTH:	465 nm
ANGLE:	100°
SUPPLY VOLTAGE:	100÷240V
BATTERY:	4.8V 2000mAh
CURRENT FLUCTUATIONS:	+/- 20%
FREQUENCY:	50÷60 Hz
PROTECTION AGAINST ELECTRIC SHOCK:	CLASS 2
INTERMITTENT OPERATION :	max 15 min
CONDITIONS OF USE:	Temperature +18°C /+40°C (+64°F /+ 140°F) Relative Humidity: < 80%
PROTECTION	IP 20
CLASS MEDICAL DEVICE:	Classe I Rule 1 93/42/CEE
SHIPPING AND STORAGE CONDITIONS:	+5/+65 C°(+41/+149 °F); 20-95% RH NON CONDENSING; 500-1060 mBar

BY DENTAL ® s.r.l.

Via Vecchia Prov. Lucchese, 49/FG • 51030 Serravalle Pistoiese (PT) - Italy • CC.I.A.A. PT n° 155715
Tel (+39)0573994355 • Fax (+39)0573919480 • CF/P.IVA 01508680475

www.bydental.com - info@bydental.com