




## ***Aspirator tips***

by **ROEKO**

# Aspire to perfection



**STERILE**



## ***Surgitip-micro and Surgitip***

Sterile, surgical aspirator tips from ROEKO, the right tip for every kind of surgical treatment.

Precise and effective aspiration as well as absolute sterility are a must for surgical procedures in the dental practice.

*Surgitip* and *Surgitip-micro* have been specially designed for use in micro surgery, endodontics and general oral surgery and fulfil exactly these requirements.

Ideal shape and length:

- The extra length of the aspirator tips means they can be held at a distance from the operating area
- The fine tube and the extra fine green tip ensure effective aspiration even in small spaces
- The angle of the tube and the conicity of the green tip ensure that the field of vision is not blocked
- The green tip provides an excellent contrast to the operation area
- As *Surgitip* and *Surgitip-micro* are made of plastic they can be used in electrosurgery

*Surgitip-micro* is specially designed for use with microscopes and loupes. With an extra fine opening of 1.2 mm the sterile tips are also ideal for use in endodontics.

*Surgitip* is available with 2.5 mm and 4.8 mm openings and are ideal for use in standard surgical procedures where sterility and effective aspiration are essential.

Ready for use in sterile peel-packs the tips are hygienic and time-saving.



Each box of *Surgitip-micro* and *Surgitip* contains an autoclavable Double Adapter for adapting the tips ( $\varnothing 6,5$ mm) to  $\varnothing 11$  mm and  $\varnothing 16$  mm aspirating equipment.



### **Aspirator tips**

These plastic *aspirator tips* with an opening of 4.8 mm are ideal for moisture evacuation during all standard procedures, where sterility is not a requirement. The length and shape correspond to those of Surgitip 4.8 mm and each pack of 100 pieces contains an autoclavable Double Adapter.

### **Adapter + Double Adapter**

To ensure that ROEKO's aspirator tips can be used on all standard aspirating equipment adapters are available. The *Double Adapter* can be used to fit *Surgitip-micro*, *Surgitip* and *Aspirator tips* onto 11 mm or 16 mm equipment (6.5 > 11 > 16 mm) while the individual Adapter can be used to fit *Concorde II* 11mm onto 16 mm equipment.



### **Concorde II**

The smooth contours of *Concorde II* aspirator tips mean they are comfortable for the patient and easy to clean. *Concorde II* is made from high quality, durable plastic.

The *Concorde II* tips ensure effective aspiration with a low noise level. The rounded flange gently retracts the tongue and cheek without causing pain to the patient. The soft contours and the smooth inner surface prevent impurities from collecting on the tips, which can be sterilised in the autoclave or chemiclave up to 134 °C. The light grey, matt surface prevents light reflection.




*Adapter* and *Double Adapter* can be sterilised in the autoclave or chemiclave up to 134 °C.

*Concorde II* is supplied in transparent boxes of 10 pieces. The tips are available with 11 mm and 16 mm diameters for all standard equipment. The 11 mm tips are supplied with an autoclavable Adapter (11 mm > 16 mm).





# roeko for effective moisture control



**Flexi Dam non latex**  
extremely elastic Dental Dam contracts closely around the tooth for best isolation. Dental Dam accessories for easy placement

**Cotton pellets**  
highly absorbent, soft and uniform



**Luna dental rolls**  
absorb moisture - throughout treatment while keeping their shape



**Comprecap anatomic**  
anatomically shaped compression caps for a dry sulcus

## Order form

		Tip opening ø mm	End opening ø mm	Contents	REF	Quantity
<b>Surgitip-micro</b>	Standard pack	1,2	6,5	20 *	462 012	
<b>Surgitip</b>	Standard pack	2,5	6,5	20 *	462 025	
	Clinic pack			100 **	462 026	
	Standard pack	4,8	6,5	20 *	462 048	
	Clinic pack			100 **	462 049	
<b>Aspirator tips</b>		4,8	6,5	100 *	461 048	
<b>Concorde II</b>			11	10 ***	450 211	
			16	10	450 216	
<b>Adapter</b>			11	1	463 025	
			16	1	463 011	
<b>Double Adapter</b>			11	1 Set	463 000	
			16			

\*incl. 1 Double Adapter    \*\*incl. 5 Double Adapters    \*\*\*incl. 1 Adapter

www.coltenewhaledent.com  
© 2004 Coltène/Whaledent



**Coltène/Whaledent Ltd.**  
President Suite Kendal House  
Victoria Way Burgess Hill  
West Sussex, RH15 9NF/U.K.  
Tel. +44 (0)1444 235486  
Fax +44 (0)1444 870640  
hq@coltenewhaledent.co.uk

**Coltène/Whaledent GmbH + Co. KG**  
Raiffeisenstrasse 30  
89129 Langenau/Germany  
Tel. +49 (0)7345 805 0  
Fax +49 (0)7345 805 201  
productinfo@coltenewhaledent.de

**Coltène/Whaledent AG**  
Feldwiesenstrasse 20  
9450 Altstätten/Switzerland  
Tel. +41 (0)71 757 53 00  
Fax +41 (0)71 757 53 01  
office@coltenewhaledent.ch

**Coltène/Whaledent Inc.**  
235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223/USA  
Tel. +1 330 916 8800  
Fax +1 330 916 7077  
productinfo@coltenewhaledent.com

