

roeko

Flexi Dam non latex

Highly elastic
Latex free
Excellent contrast



roeko

Flexi Dam non latex

Highly elastic



Flexi Dam non latex is an exceptionally elastic dental dam. It is made from an elastic plastomer and can be elongated more than 1000 % before tearing. This means it is more tenacious than latex dam and is really simple to place. Flexi Dam non latex should be stretched before use. The punched hole should be 1 or 2 sizes smaller than normally used for latex dam.

Latex free



As the dam is latex free and powder free, patients and dental staff do not come into contact with allergens when Flexi Dam non latex is used. It contracts extremely well around the cervix of tooth for effective isolation and an absolutely dry working area.

Excellent contrast

Flexi Dam non latex is available in blue and violet. The violet colour is unique for dental dam and can easily be distinguished from latex materials.



The blue Flexi Dam non latex provides a good contrast to the working area and has a calming effect on the user. These are important advantages when performing aesthetic dentistry and extensive treatment.

Flexi Dam non latex

- more tear resistant than latex
- highly elastic
- excellent recovery properties
- powderfree
- latex free
- no unpleasant smell


Material: elastic plastomer
Colour: purple, blue
Measurements: approx. 150 x 150 mm
Thickness: approx. 0.50 mm

Order information:

Flexi Dam non latex	30 pcs. + template
purple	REF 390 035
blue	REF 390 033

Roeko Dental Dam Set	REF 390 000
Punch Type Ainsworth	1 pc.
Clamp forceps Type Ivory	1 pc.
Clamp Set	6 clamps
Flexi Dam non latex purple	30 pcs. + template
Metal frame acc. to Young	1 pc.

www.coltenewhaledent.com
© 2006 Coltène/Whaledent AG



Coltène/Whaledent Ltd.
 The President Suite-A
 Kendal House, Burgess Hill
 Victoria Way
 West Sussex, RH15 9NF/U.K.
 Tel. +44 (0)1444 235486
 Fax +44 (0)1444 870640
 hq@coltenewhaledent.co.uk

Coltène/Whaledent GmbH + Co. KG
 Raiffeisenstrasse 30
 89129 Langenau/Germany
 Tel. +49 (0)7345 805 0
 Fax +49 (0)7345 805 201
 productinfo@coltenewhaledent.de

Coltène/Whaledent Inc.
 235 Ascot Parkway
 Cuyahoga Falls, Ohio 44223/USA
 Tel. +1 330 916 8800
 Fax +1 330 916 7077
 productinfo@coltenewhaledent.com

Coltène/Whaledent AG
 Feldwiesenstrasse 20
 9450 Altstätten/Switzerland
 Tel. +41 (0)71 757 53 00
 Fax +41 (0)71 757 53 01
 office@coltenewhaledent.ch

coltène whaledent