

Detartrine

- Last Updated (Thursday, 09 April 2009)

This product does not need a Safety Data Sheet/COSHH.

Detartrine

Abrasive paste for scaling

Formula

Silica 55.00 g

35 % formaldehyde solution 0.13 g

Excipient q.s.ad 100 g Properties

Détartrine is presented as pink paste of very pleasant taste.

Détartrine has a purely mechanical abrasive action able to perfect
ultra-sounds, without eroding the dental enamel.

The essential oils present in Détartrine produce a tingling fresh
patient. Indications

the scaling, after having used scalers or

sensation which is greatly appreciated by the

- Finishing of scaling
- Elimination of colourations due to tobacco or food
- Fine polishing of fillings.

Directions for use

Using a spatula, place a sufficient quantity of Détartrine on the
manner. Precautions for use

Possible risk of allergic reactions to formaldehyde. Remark

Rinsing may be carried out during or at the end of the treatment.

Jar containing 45 g.

circular brush, and proceed in the usual

Presentations

For professional use only in the practice of dentistry.