

Raypex® 5

Apex Locator – for successful endodontic treatments



5 good reasons!




Accurate length determination – the foundation for successful endodontic treatments

→ Reliable 1

- Large projection of the apical section with unique visualization of the apical constriction
- Pre-calibrated display of the apical constriction area, no individual calibration necessary¹
- Live representation of the file movement along the entire root canal

Reliable
and accurate



→ Accurate 2

- Impedance measurement based on advanced multi-frequency system
- Latest digital technology





→ User-friendly **3**

- Foldable display adjusts to desired viewing angle
- Anti-slip casing with cable storage compartment
- Rechargeable AAA batteries
- Battery charger indicator
- Instruction CD-ROM included in kit

→ Patient-friendly **4**

- Built-in demo mode for quick demonstration to patients
- Length measurement without radiation

Patient-friendly
and safe

→ Safe **5**

- Electronic length determination – safest method
No more over-extended or under-extended root fillings
- Conventional radiography diagnostics – not reliable due to the 2D representation of a 3D root canal morphology

User-friendly



100% confident?



but 100% off target!

Raypex[®] 5 for successful endodontic treatments!

¹ study data on file



Raypex® 5



Raypex® 5 Kit

Ref. V04 0145 000 511 Raypex® 5 Kit

1 Raypex® 5 Apex locator
1 measuring cable, 1 charger
2 lip clips, 2 file clips, 1 touch probe
1 instruction CD-ROM



Accessories

Ref. V04 0141 000 506 1 measuring cable

Ref. V04 0141 000 507 1 battery

Ref. V04 0141 000 514 5 lip clips

Ref. V04 0141 000 515 2 file clips

Ref. V04 0141 000 518 2 touch probes

VDW GmbH • P.O. Box 830954 • 81709 Munich, Germany
Tel. +49 89/627 34-0 • Fax +49 89/627 34-195
info@vdw-dental.com • www.vdw-dental.com

VW000188 Rev. 3/11.8.05

**VDW**[®]
ENDODONTIC SYNERGY

**Endo
Easy
Efficient**[®]

