

Admira® Flow

Flowable light-curing



Ormocer®-based restorative

Indications

Filling of small cavities

Extended fissure sealing

Blocking-out of undercuts

Lining, coating of cavities (lining technique, CBF-technique)

Cavities class III and V

Repair of composite fillings and veneers

Luting of full ceramic restorative pieces (porcelain-only-crowns, porcelain inlays)

Advantages

- Optimal flowability (thixotropic): fine-flowing during application, stable in the state of rest
- Stress-breaking elasticity
- Exact match of translucency and shade to the packable Admira
- Fine-flowing with excellent wetting capability
- In syringes or Caps with small long cannulae
- Biocompatibility
- Low shrinkage
- Time-saving application

Presentation

Art.-No. 2481 Triset 3 × 1,8 g syringes (A2, A3, A3.5), 5 ml Vococid gel, 2 ml Admira Bond, accessories

Refills 2 × 1,8 g syringes, shades:

Art.-No. 2486	A1
Art.-No. 2482	A2
Art.-No. 2483	A3
Art.-No. 2484	A3.5
Art.-No. 2487	A4
Art.-No. 2488	OA3.5