

Futurabond NR

SingleDose

Light-curing nano-reinforced self-etch dentine/enamel bond



Indication

Direct restoration of all cavity classes with all light-curing restoratives (especially nanocomposites)

Treatment of dentine hypersensitivity

Advantages

- Excellent adhesion to dentine and enamel
- Durable marginal integrity
- Moisture-tolerant
- Easy and fast application - one-coat - one-cure technique
- Universal for all light-curing restoratives
- Ideal for paediatric use
- Fluoride release
- No refrigerated storage needed



k and hygienic - one-step activation blister

- No additional devices required
- No spilling in any position

Presentation

Art.-Nr. 1152 50 × *SingleDose*, accessories

Art.-Nr. 1153 200 × *SingleDose*, accessories