

# Grandio® Flow

## Flowable universal nano-hybrid restorative



### Indications

- Filling of smaller cavities
- Extended fissure sealing
- Blocking-out of undercuts
- Lining, coating of cavities (CBF technique, lining technique)
- Restorations class III to V
- Repair of fillings and veneers
- Luting of translucent prosthetic pieces (e.g. porcelain-only crowns, veneers)

### Advantages

- Optimal flowability and high filler degree
- Exact shade match to Grandio
- Excellent wetting properties
- Outstanding physical properties
- Significantly lower polymerization shrinkage than conventional flowables
- Very good abrasion resistance
- High transverse strength
- Can be used with all conventional bonding agents
- The shades BL (bleach light) and WO (white opaque) fulfil special indications: "Bleach light" is suited for bleached teeth and for use in paediatric dentistry. The shade "white opaque" makes an ideal restoration base on e.g. discoloured dentine areas and in core build-up.

### Presentation

Refill packs 2 × 1.8 g syringes, shades:

Art.-No. 1863 A1

Art.-No. 1864 A2  
Art.-No. 1865 A3  
Art.-No. 1866 A3.5  
Art.-No. 1868 A4  
Art.-No. 1867 B2  
Art.-No. 1873 C2  
Art.-No. 1874 D2  
Art.-No. 1869 B1  
Art.-No. 1870 OA2  
Art.-No. 1871 BL (bleach light)  
Art.-No. 1872 WO (white opaque)